EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Any serious incident relating to the use of any device should be reported to the manufacturer and Competent Authority of the Member State in which the user or patient is located.

Application Guidelines

- Use a small test patch if person has known allergy to other adhesives
- Apply to clean, dry skin at least 60 minutes before exercise, swimming or showering
- Remove hair
- Round the corners of the tape to reduce lifting or peeling
- Tear paper to peel backing paper to avoid handling the adhesive
- Apply one end of the tape (anchor point) with no tension. This should be at least 2" long and longer if more force is expected (e.g. PowerBand, fast movements or lots of stretch)
- Position the body part in the shortened position
- Hold the anchor point and take up the slack in the tape (do not stretch strongly).
- Maintain pressure on the anchor as you stick and smooth the tape down
- Apply the other end (anchor) with no tension
- PRESS & HOLD the ends for 30 seconds to 1 minute and then rub thoroughly to activate the glue
- Avoid starting or finishing one piece of tape on another and limit overlapping due to the strong shearing forces.

Tape Removal

- Hold the skin adjacent to the tape and peel the tape back along itself in short stages. Try to pull the skin away from the tape at the same time.
- If irritation has occurred a coconut oil or similar may be used to aid removal.

Warning

- Remove immediately if any itching, burning, stinging, irritation, pins and needles or numbness or a deterioration in symptoms occurs.
- Always consult a healthcare professional for advice
- This information is for education purposes only and does not constitute medical advice for any particular individual or case.

Caution

- Do not use on frail, damaged, open or sunburned skin or where there is circulatory compromise
- Traction blisters may result from shearing on skin due to incorrect application
- Although uncommon, allergic reactions may occur

CAUTION

This product requires careful and skilful application. Failure to follow the Directions for Use available on the product packaging or at www.dynamictape.com can result in skin irritation, blisters and poor adhesion.

1 - Remove hair

2 - Clean & Dry (e.g. alcowipe)


3 - Rub to heat skin


4 - Apply Adhesive Spray*

5 - Anchor Point – No Stretch

6 - Anchor with thumb and tension in opposite direction to tape. This will remove skin folds and reduce tension on the skin at the end of the tape

7- Apply Gentle Tension

8- Apply other anchor point with no tension

* Very good results can be obtained without the use of an adhesive spray provided that all other guidelines are adhered to and sufficient time (> one hour) is allowed before participating in vigorous exercise, swimming or bathing or strongly stretching the tape.

9 - Rub thoroughly to activate glue,
PRESS & HOLD the ends for 30-60 seconds


10 - Allow 60-90 minutes before engaging in
vigorous exercise, swimming or bathing


11 - Always use spray on feet and ankles*


12 - Spray back of first layer before applying the
next layer when overlapping tape*

13 - Always lock over foot and achilles
to increase load absorption as well as adhesion


14 - Cover and protect

WARN – Remove tape immediately if itching, stinging, burning or irritation occurs as you may be developing a reaction which can lead to skin breakdown.

Poor application will result in poor adhesion, tension, shearing and blisters

(c) 2012 PosturePals Pty Ltd
All rights reserved. No part of this publication may be reproduced, stored in any retrieval system or transmitted, in any form or by any means, electronic, mechanical, photocopying, recording or otherwise without the prior permission of the publishers.

CAUTION

This product requires careful and skilful application. Failure to follow the Directions for Use available on the product packaging or at www.dynamictape.com can result in skin irritation, blisters and poor adhesion.

1 - Remove hair

2 - Clean & Dry (e.g. alcowipe)


3 - Rub to heat skin


4 - Apply Adhesive Spray*

5 - Anchor Point – No Stretch

6 - Anchor with thumb and tension in opposite direction to tape. This will remove skin folds and reduce tension on the skin at the end of the tape

7- Apply Gentle Tension

8- Apply other anchor point with no tension

* Very good results can be obtained without the use of an adhesive spray provided that all other guidelines are adhered to and sufficient time (> one hour) is allowed before participating in vigorous exercise, swimming or bathing or strongly stretching the tape.

9 - Rub thoroughly to activate glue,
PRESS & HOLD the ends for 30-60 seconds


10 - Allow 45-60 minutes before engaging in vigorous exercise, swimming or bathing


11 - Always use spray on feet and ankles*


12 - Spray back of first layer before applying the next layer when overlapping tape*

13 - Always lock over foot and achilles to increase load absorption as well as adhesion

14 - Cover and protect

WARN – Remove tape immediately if itching, stinging, burning or irritation occurs as you may be developing a reaction which can lead to skin breakdown.

Poor application will result in poor adhesion, tension, shearing and blisters

(c) 2012 PosturePals Pty Ltd
All rights reserved. No part of this publication may be reproduced, stored in any retrieval system or transmitted, in any form or by any means, electronic, mechanical, photocopying, recording or otherwise without the prior permission of the publishers.

Contra-indications

- Do not use on frail, broken or sunburned skin or on those with known sensitivity

Caution

- Apply strictly according to directions for use
- Accredited Dynamic Tape training recommended prior to use
- DO NOT stretch strongly. Excessive stretch will result in blisters
- Remove immediately, discontinue use and seek medical advice if itching, burning, stinging, rash, redness or irritation occurs
- Circumferential applications should be applied on an angle to prevent compression of blood vessels and nerves
 - Remove if pins and needles or numbness occurs
- Check product thoroughly if packaging is damaged


Application

- Remove hair (clippers recommended)
- Clean and dry skin (remove creams, lotions & oils)
- Round off the corners of the tape to reduce lifting
- Place the body part in desired position
- Apply an anchor point of > 4cm with no stretch to avoid tension on skin
- Hold anchor (4cm from end) to minimise traction on the skin and gently stretch tape until resistance is first felt – DO NOT stretch strongly
- Apply the final 4cm of tape with no stretch
- Rub thoroughly to activate the heat sensitive glue
- Allow 45 – 60 minutes before swimming, showering or vigorous exercise
- May remain in place for up to five days as directed
- Remove tape in the direction of hair growth. Hold down skin and peel tape back along itself
- Do not remove when wet


Further Tips

- Use adhesive spray to hasten bonding, around foot and ankle or if overlapping tape (apply on back of first layer)
- Apply a locking strip over foot and ankle applications to improve adhesion and load absorption
- A rigid locking strip (zinc oxide) may be useful around the ends of the tape. This is especially useful on fingers


For videos, specific techniques, tips and disclaimer please visit www.dynamictape.com

A


B


C

PosturePals Box applications:


- A - Thoracic Box (vertical application)
- B - Thoracic Box (horizontal application)
- C - Lumbar Box


A


B

C


D

E

PosturePals X applications:

- A - Pregnancy X
- B - Knee X
- C - Lumbar X
- D - Inhibitory X
- E - Thoracic X
- F - Shoulder X

F

English - En

Application Guidelines

1. Apply in the position as shown.
2. Ensure that the skin is clean and dry and free from hair, creams and oils.
3. Do not use if you have sunburn, frail or broken skin or known sensitivity to sticking plasters or strapping tape.
4. Slowly and carefully remove one section of the backing paper.
5. Apply one quarter of the PosturePals tape at a time before removing the next section of backing paper.
6. Stretch the tape until some resistance is felt. DO NOT stretch strongly
7. Ensure that the ends of the tape are not stretched and remain circular to reduce the risk of blisters caused by the tape pulling on the skin.
8. Remove immediately and seek medical advice if necessary if any itching, burning, stinging or irritation occurs.
9. PosturePals can generally remain in place for four to five days and will tolerate brief showers.
10. Remove by holding the skin down and gently peeling the tape back on itself.

If you have spinal stenosis, moderate to severe scoliosis or are unsure if this device is appropriate for you, please consult with your health care practitioner prior to use.

Please visit www.posturepals.com prior to use for instructional videos and further information including our waiver/disclaimer and answers to frequently asked questions.

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