English - En

Application Guidelines

1. Apply in the position as shown.
2. Ensure that the skin is clean and dry and free from hair, creams and oils.
3. Do not use if you have sunburn, frail or broken skin or known sensitivity to sticking plasters or strapping tape.
4. Slowly and carefully remove one section of the backing paper.
5. Apply one quarter of the PosturePals tape at a time before removing the next section of backing paper.
6. Stretch the tape until some resistance is felt. DO NOT stretch strongly
7. Ensure that the ends of the tape are not stretched and remain circular to reduce the risk of blisters caused by the tape pulling on the skin.
8. Remove immediately and seek medical advice if necessary if any itching, burning, stinging or irritation occurs.
9. PosturePals can generally remain in place for four to five days and will tolerate brief showers.
10. Remove by holding the skin down and gently peeling the tape back on itself.

If you have spinal stenosis, moderate to severe scoliosis or are unsure if this device is appropriate for you, please consult with your health care practitioner prior to use.

Please visit www.posturepals.com prior to use for instructional videos and further information including our waiver/disclaimer and answers to frequently asked questions.

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